

Goat Anti-Mouse IgG(H+L)(Elab Fluor[®] 488 conjugated)

Catalog Number:E-AB-1056

6 Publications



Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Reactivity	Mouse
Immunogen	Mouse IgG
Host	Goat
Purification	Affinity purification
Conjugation	Elab Fluor [®] 488 (Ex/Em:493nm/519nm)
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4

Applications Recommended Dilution

Immunofluorescence	1:100-300
---------------------------	-----------

Preparation & Storage

Storage	Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.
----------------	---

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017