

## Influenza A NP (366-374)Strain A/NT/60/68

**Catalog No.** E-PP-1380

**Storage** Shipped at 4°C. Stored at -20°C for one year.

**Applications**

**Reactivity**

**Host**

**Isotype**

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.

### Images

### Immunogen Information

### Product Information

#### For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit [www.elabscience.com](http://www.elabscience.com).

**Focus on your research  
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.