

Melanocyte Associated Antigen gp 100 (17-25)

Catalog No. E-PP-1529

Storage Shipped at 4°C. Stored at -20°C for one year.

Applications

Reactivity

Host

Isotype

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images

Immunogen Information

Product Information

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.