

Platelet Factor 4 (58-70), human

| | | |
|---------------------|---|-------------------|
| Catalog No. | E-PP-1773 | Reactivity |
| Storage | Shipped at 4°C. Stored at -20°C for one year. | Host |
| Applications | | Isotype |

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images

Immunogen Information

Product Information

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.