

Rabbit Anti-Goat IgG(H+L)(CY3 conjugated)

Catalog No. E-AB-1012

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Product Information

Immunogen: Goat IgG

Purification Method: Affinity purification

Conjugation: CY3

Concentration: 0.2mg/mL(specific antibody purification≥99%),CY3:Ab=3-9:1

Applications: IF

Dilution: Immunofluorescence: 1:50-100

Optimal working dilutions must be determined by end user. Dilute before when it's needed, and the residual should be discarded.

Amax: 554nm

E_{max}: 568-574nm

Buffer: PBS with 0.05% sodium azide and 50% glycerol, PH7.4

Storage: Store at -20°C. Avoid freeze / thaw cycles. Protected from prolonged exposure to light.