

Recombinant Human IL4RA/CD124 Protein (mFc Tag)

Catalog No. PKSH033612

Description

Synonyms	Interleukin-4 receptor subunit alpha, IL-4 receptor subunit alpha, IL-4R subunit alpha, IL-4R-alpha, IL-4RA, CD124, IL-4-binding protein, IL4-BP, IL4R, IL4RA
Species	Human
Expression_host	HEK293 Cells
Sequence	Met26-Gln231
Accession	P24394
Mol_Mass	50.2 kDa
AP_Mol_Mass	65-86 kDa
Tag	C-mFc
Bio_activity	Testing in progress

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS:pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Background

Interleukin 4 Receptor alpha (IL4-Ra) is a widely expressed 140 kDa transmembrane glycoprotein in the class I cytokine receptor family. Mature human IL4-Ra consists of a 207 amino acid (aa) extracellular domain (ECD) that contains a cytokine binding region and one fibronectin type III domain; a 24 aa transmembrane segment; and a 569 aa cytoplasmic domain that contains one Box 1 motif and one ITIM motif. IL4-Ra plays an important role in Th2-biased immune responses; alternative macrophage activation; mucosal immunity; allergic inflammation; tumor progression; and atherogenesis. Soluble forms of IL4-Ra; generated by alternate splicing or proteolysis; retain ligand binding properties and inhibit IL-4 bioactivity.

SDS-PAGE

