

Recombinant Human Interleukin-2/IL-2 Protein (C145S)(Active)

Catalog No. PKSH033623

Description

Synonyms	Interleukin-2, IL-2, T-Cell Growth Factor, TCGF, Aldesleukin, IL6
Species	Human
Expression_host	E.coli
Sequence	Pro22-Thr153(Cys145Ser)
Accession	P60568
Mol_Mass	15.5 kDa
AP_Mol_Mass	14 kDa
Tag	None
Bio_activity	Measured in a cell proliferation assay using CTLL-2 mouse cytotoxic T cells. The ED50 for this effect is 90-270 pg/ml

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 10mM Acetata-Na, 5% Trehalose, pH 4.5. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the p
Reconstitution	Please refer to the printed manual for detailed information.

Background

Interleukin-2(IL-2) is an interleukin; a type of cytokine signaling molecule in the immune system;belongs to the IL-2 family. It is a powerful immunoregulatory lymphokine produced by T-cells in response to antigenic or mitogenic stimulation. IL-2/IL-2R signaling is required for T-cell proliferation and other fundamental functions that are essential for the immune response. IL-2 stimulates growth and differentiation of B-cells; NK cells; lymphokine-activated killer cells; monocytes; macrophages and oligodendrocytes.

SDS-PAGE

