


TSHB Polyclonal Antibody


Catalog No.	E-AB-40009	Reactivity	M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB,IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Images



Western Blot analysis of Rat lung and Mouse lung tissue using TSHB Polyclonal Antibody at dilution of 1:600



Immunohistochemistry of paraffin-embedded Rat lung using TSHB Polyclonal Antibody at dilution of 1:50



Immunohistochemistry of paraffin-embedded Rat brain using TSHB Polyclonal Antibody at dilution of 1:50

Immunogen Information

Immunogen	Recombinant Rat Thyrotropin subunit beta protein
Swissprot	P04652
Synonyms	CHNG4, TSHB

Product Information

Calculated MW	16kDa
Buffer	PBS with 0.05% Proclin300 and 50% glycerol, pH7.4.
Purify	Affinity purification
Dilution	WB 1:500-1:1000, IHC 1:50-1:100

Background

Indispensable for the control of thyroid structure and metabolism.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.