


FABP3 Polyclonal Antibody

Catalog No.	E-AB-40085	Reactivity	R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Rat heart using FABP3 Polyclonal Antibody at dilution of 1:50

Immunogen Information

Immunogen	Recombinant Rat Fatty acid-binding protein, heart protein
Swissprot	P07483
Synonyms	FABP 11,FABP 3,FABP11,FABP3,FABPH,Fatty acid binding protein 3,Fatty acid binding protein 3 muscle and heart H-FABP,M FABP,M-FABP,MDGI

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IHC 1:50-1:100

Background

FABP are thought to play a role in the intracellular transport of long-chain fatty acids and their acyl-CoA esters.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.