


TNP2 Polyclonal Antibody

Catalog No.	E-AB-61354	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of U2OS cells using TNP2 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human TNP2 (NP_005416.1).
GeneID	7142
Swissprot	Q05952
Synonyms	TNP2, TP2

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

lays a key role in the replacement of histones to protamine in the elongating spermatids of mammals. In condensing spermatids, loaded onto the nucleosomes, where it promotes the recruitment and processing of protamines, which are responsible for histone eviction.
STP2_HUMAN,Q05952

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.