


PHYHD1 Polyclonal Antibody

Catalog No.	E-AB-61585	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of A-549 cells using PHYHD1 Polyclonal Antibody

Immunogen Information

Immunogen	Recombinant fusion protein of human PHYHD1 (NP_001094346.1).
GeneID	254295
Swissprot	Q5SRE7
Synonyms	PHYHD1

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

PHYHD1 (Phytanoyl-CoA Dioxygenase Domain Containing 1) is a Protein Coding gene. Gene Ontology (GO) annotations related to this gene include dioxygenase activity.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.