


DAAM2 Polyclonal Antibody

Catalog No.	E-AB-61723	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Rat brain using DAAM2 Polyclonal Antibody at dilution of 1:100 (40x lens).

Immunogen Information

Immunogen	Recombinant fusion protein of human DAAM2 (NP_056160.2).
GeneID	23500
Swissprot	Q86T65
Synonyms	DAAM2,dJ90A20A.1

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IHC 1:50-1:200

Background

DAAM2 (Dishevelled Associated Activator Of Morphogenesis 2) is a Protein Coding gene. Diseases associated with DAAM2 include Molybdenum Cofactor Deficiency and Molybdenum Cofactor Deficiency, Complementation Group A. Among its related pathways are Wnt signaling pathway and WNT Signaling. Gene Ontology (GO) annotations related to this gene include binding and Rho GTPase binding. An important paralog of this gene is DAAM1.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.