


## SPAG16 Polyclonal Antibody


<b>Catalog No.</b>	E-AB-61737	<b>Reactivity</b>	H,M,R
<b>Storage</b>	Store at -20°C. Avoid freeze / thaw cycles.	<b>Host</b>	Rabbit
<b>Applications</b>	WB,IHC	<b>Isotype</b>	IgG

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


### Images



Western blot analysis of extracts of Mouse testis using SPAG16 Polyclonal Antibody at dilution of 1:1000.



Immunohistochemistry of paraffin-embedded Human brain cancer using SPAG16 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Mouse brain using SPAG16 Polyclonal Antibody at dilution of 1:100 (40x lens).

### Immunogen Information

<b>Immunogen</b>	Recombinant fusion protein of human SPAG16 (NP_001020607.1).
<b>GeneID</b>	79582
<b>Swissprot</b>	Q8N0X2
<b>Synonyms</b>	SPAG16,PF20,WDR29

### Product Information

<b>Calculated MW</b>	20kDa/29kDa/39kDa/50kDa/70kDa
<b>Observed MW</b>	71kDa
<b>Buffer</b>	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
<b>Purify</b>	Affinity purification
<b>Dilution</b>	WB 1:500-1:2000 IHC 1:50-1:200

### Background

SPAG16 (Sperm Associated Antigen 16) is a Protein Coding gene. Diseases associated with SPAG16 include Pulmonary Subvalvular Stenosis and Primary Ciliary Dyskinesia. Gene Ontology (GO) annotations related to this gene include protein kinase binding.

#### For Research Use Only

Thank you for your recent purchase.  
 If you would like to learn more about antibodies, please visit [www.elabscience.com](http://www.elabscience.com).

**Focus on your research  
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.