


SPAG16 Polyclonal Antibody


Catalog No.	E-AB-61737	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB,IHC	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Images



Western blot analysis of extracts of Mouse testis using SPAG16 Polyclonal Antibody at dilution of 1:1000.



Immunohistochemistry of paraffin-embedded Human brain cancer using SPAG16 Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunohistochemistry of paraffin-embedded Mouse brain using SPAG16 Polyclonal Antibody at dilution of 1:100 (40x lens).

Immunogen Information

Immunogen	Recombinant fusion protein of human SPAG16 (NP_001020607.1).
GeneID	79582
Swissprot	Q8N0X2
Synonyms	SPAG16,PF20,WDR29

Product Information

Calculated MW	20kDa/29kDa/39kDa/50kDa/70kDa
Observed MW	71kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000 IHC 1:50-1:200

Background

SPAG16 (Sperm Associated Antigen 16) is a Protein Coding gene. Diseases associated with SPAG16 include Pulmonary Subvalvular Stenosis and Primary Ciliary Dyskinesia. Gene Ontology (GO) annotations related to this gene include protein kinase binding.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.