


GNG7 Polyclonal Antibody

Catalog No.	E-AB-62770	Reactivity	M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of various cell lines using GNG7 Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	Recombinant fusion protein of human GNG7 (NP_443079.1).
GeneID	2788
Swissprot	O60262
Synonyms	GNG7

Product Information

Calculated MW	7kDa
Observed MW	12kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:1000

Background

GNG7 (G Protein Subunit Gamma 7) is a Protein Coding gene. Diseases associated with GNG7 include Keratinizing Squamous Cell Carcinoma and Esophageal Cancer. Among its related pathways are Translation Translation regulation by Alpha-1 adrenergic receptors and ADP signalling through P2Y purinoceptor 12. Gene Ontology (GO) annotations related to this gene include obsolete signal transducer activity. An important paralog of this gene is GNG12.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.