


ZSCAN4C Polyclonal Antibody

Catalog No.	E-AB-63298	Reactivity	M
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western blot analysis of extracts of Mouse kidney using ZSCAN4C Polyclonal Antibody at dilution of 1:1000.

Immunogen Information

Immunogen	A synthetic peptide of mouse ZSCAN4
GeneID	245109
Swissprot	Q80VJ6
Synonyms	ZSCAN4,ZNF494,ZSCAN4C

Product Information

Calculated MW	57kDa
Observed MW	57kDa
Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	WB 1:500-1:2000

Background

ZSCAN4 (Zinc Finger And SCAN Domain Containing 4) is a Protein Coding gene. Diseases associated with ZSCAN4 include Facioscapulohumeral Muscular Dystrophy 1 and Facioscapulohumeral Muscular Dystrophy 2. Among its related pathways are Preimplantation Embryo. Gene Ontology (GO) annotations related to this gene include DNA-binding transcription factor activity. An important paralog of this gene is ZNF263.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.