


RHBDD1 Polyclonal Antibody

Catalog No.	E-AB-64986	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunofluorescence analysis of L929 cells using RHBDD1 Polyclonal Antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen	Recombinant fusion protein of human RHBDD1 (NP_115652.2).
GeneID	84236
Swissprot	Q8TEB9
Synonyms	RHBDD1,RRP4

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IF 1:50-1:100

Background

RHBDD1 (Rhomboid Domain Containing 1) is a Protein Coding gene. Diseases associated with RHBDD1 include Benign Mammary Dysplasia and Brain Glioblastoma Multiforme. Gene Ontology (GO) annotations related to this gene include serine-type endopeptidase activity and endopeptidase activity. An important paralog of this gene is RHBDD3.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.