


HIGD1A Polyclonal Antibody


Catalog No.	E-AB-65044	Reactivity	H
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	IHC,IF	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Immunohistochemistry of paraffin-embedded Human colon carcinoma using HIGD1A Polyclonal Antibody at dilution of 1:100 (40x lens).



Immunofluorescence analysis of U-2 OS cells using HIGD1A Polyclonal Antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.

Immunogen Information

Immunogen	A synthetic peptide of human HIGD1A (NP_054775.2).
GeneID	25994
Swissprot	Q9Y241
Synonyms	HIGD1A,HIG1,RCF1a

Product Information

Buffer	PBS with 0.02% sodium azide, 50% glycerol, pH7.3.
Purify	Affinity purification
Dilution	IHC 1:50-1:100 IF 1:50-1:100

Background

HIGD1A (HIG1 Hypoxia Inducible Domain Family Member 1A) is a Protein Coding gene. Diseases associated with HIGD1A include Hypoxia and Behr Syndrome. Among its related pathways are Regulation of activated PAK-2p34 by proteasome mediated degradation and Cellular Senescence (REACTOME). An important paralog of this gene is HIGD1C.

For Research Use Only

Thank you for your recent purchase.
 If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.