


Phospho-PDGFRα (Tyr849) Polyclonal Antibody

Catalog No.	E-AB-21070	Reactivity	H,M,R
Storage	Store at -20°C. Avoid freeze / thaw cycles.	Host	Rabbit
Applications	WB,IHC-p,ELISA	Isotype	IgG

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Images



Western Blot analysis of 293 cells with Phospho-PDGFRα (Tyr849) Polyclonal Antibody

Immunogen Information

Immunogen	Synthesized peptide derived from human PDGFR-α around the phosphorylation site of Y849.
Swissprot	P16234
Synonyms	PDGFRA,PDGFR2,RHEPDGFRA,Platelet-derived growth factor receptor alpha,PDGF-R-alpha,PDGFR-alpha,Alpha platelet-derived growth factor receptor,Alpha-type platelet-derived growth factor receptor,CD140 antigen-like family member A,CD14

Product Information

Calculated MW	123kDa
Observed MW	140kDa
Buffer	PBS with 0.02% sodium azide, 0.5% protective protein and 50% glycerol, pH7.4
Purify	Affinity purification
Dilution	WB 1:500-1:2000, IHC 1:100-1:300, ELISA 1:5000

Background

Receptor that binds both PDGFA and PDGFB and has a tyrosine-protein kinase activity.

For Research Use Only

Thank you for your recent purchase.

If you would like to learn more about antibodies, please visit www.elabscience.com.

**Focus on your research
 Service for life science**

Applications:WB-Western Blot IHC-Immunohistochemistry IF-Immunofluorescence IP-Immunoprecipitation FC-Flow cytometry ChIP-Chromatin Immunoprecipitation Reactivity: H-Human R-Rat M-Mouse Mk-Monkey Dg-Dog Ch-Chicken Hm-Hamster Rb-Rabbit Sh-Sheep Pg-Pig Z-Zebrafish X-Xenopus C-Cow.