

# Recombinant Human Troponin C/TNNC1 Protein (His Tag)



Catalog Number:PKSH033520

**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


## Description

<b>Synonyms</b>	CMD1Z;CMH13;TN-C;TNC;TNNC
<b>Species</b>	Human
<b>Expression Host</b>	E.coli
<b>Sequence</b>	Met1-Glu161
<b>Accession</b>	P63316
<b>Calculated Molecular Weight</b>	19.8 kDa
<b>Observed molecular weight</b>	17-20 kDa
<b>Tag</b>	N-His

## Properties

<b>Purity</b>	> 95 % as determined by reducing SDS-PAGE.
<b>Endotoxin</b>	< 1.0 EU per µg of the protein as determined by the LAL method.
<b>Storage</b>	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
<b>Shipping</b>	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
<b>Formulation</b>	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 100mM NaCl, 1mM DTT, 10% Glycerol, pH 8.0.
<b>Reconstitution</b>	Not Applicable

## Data



> 95 % as determined by reducing SDS-PAGE.

## Background

Troponin I, also known as TNI, is a 24 kDa component of a protein complex on striated muscle thin filaments. Troponin is the central regulatory protein of striated muscle contraction. Tn consists of three components: Tn-I which is the inhibitor of actomyosin ATPase, Tn-T which contains the binding site for tropomyosin and Tn-C. The binding of calcium to Tn-C abolishes the inhibitory action of Tn on actin filaments. Troponin I inhibits the calcium-dependent muscle contraction mediated by Troponins C and T. The expression of cardiac Troponin I (TNNI3) is restricted to cardiac muscle, while TNNI1 and TNNI2 (encoded by distinct genes) are expressed in skeletal muscle. Mutations of cardiac Troponin I are associated with hereditary cardiomyopathy. Human cardiac Troponin I shares 93% amino acid sequence identity with mouse and rat cardiac Troponin I.

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017