

Recombinant Human Vaspin/Serpina12 Protein (GST Tag)

Catalog Number:PKSH033534

DIA • AN[®]
by Elabscience

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Serpin A12;OL-64;Visceral Adipose Tissue-Derived Serine Protease Inhibitor;Vaspin;Visceral Adipose-Specific Serpin;SERPINA12
Species	Human
Expression Host	E.coli
Sequence	Leu21-Lys414
Accession	Q8IW75
Calculated Molecular Weight	71.4 kDa
Observed molecular weight	72 kDa
Tag	N-GST

Properties

Purity	> 80 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
Formulation	Supplied as a 0.2 µm filtered solution of 50mM Tris-HCl, 160mM NaCl, 0.2mM PMSF, 1mM DTT, 10% Glycerine, pH 7.2 .
Reconstitution	Not Applicable

Data



> 80 % as determined by reducing SDS-PAGE.

Background

Vaspin (Visceral Adipose-Specific SERPIN) is a newly described adipokine. Vaspin has three β -sheets; nine α -helices; and one central loop; the structure is part of the set of distinctive features that are descriptive of Serpin family members. Vaspin is also a unique insulin sensitizing adipocytokine in obesity. A recent publication indicates that Vaspin mRNA expression in visceral fat is positively correlated with BMI and percent of body fat, and could be associated with parameters of obesity; insulin resistance; and glucose metabolism. These findings suggest a potential clinical use for Vaspin in ameliorating certain aberrations seen in the obesity metabolic syndrome.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017