

Recombinant Human RHEB Protein (GST Tag)

Catalog Number:PKSH033545



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	GTP-Binding Protein Rheb;Ras Homolog Enriched in Brain;RHEB;RHEB2
Species	Human
Expression Host	E.coli
Sequence	Met1-Met184
Accession	Q15382
Calculated Molecular Weight	20.4 kDa
Observed molecular weight	45 kDa
Tag	N-GST

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 8% Trehalose, 4% Mannitol, 50mM NaCl, 10mM GSH, 0.05% Tween 80, pH6.5. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

Background

GTP-Binding Protein Rheb (RHEB) is a member of the small GTPase superfamily and encodes a lipid-anchored; cell membrane protein with five repeats of the RAS-related GTP-binding region. Highest levels of RHEB can be found in the skeletal and cardiac muscle; and it is vital in the regulation of growth and cell cycle progression due to its role in the Insulin/TOR/S6K signaling pathway. RHEB stimulates the phosphorylation of S6K1 and EIF4EBP1 through activation of mTORC1 signaling; and it activates the protein kinase activity of mTORC1. RHEB has GTPase activity and shuttles between a GDP-bound form and a GTP-bound form; farnesylation of the protein is required for this activity.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017