

Recombinant Human MAP1LC3B Protein

Catalog Number:PKSH033588



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Microtubule-associated proteins 1A/1B light chain 3B;Autophagy-related protein LC3 B;Autophagy-related ubiquitin-like modifier LC3 B;MAP1 light chain 3-like protein 2;MAP1A/MAP1B light chain 3 B;MAP1A/MAP1B LC3 B;Microtubule-associated protein 1 light chain 3 beta;MAP1LC3B;MAP1ALC3
Species	Human
Expression Host	E.coli
Sequence	Met1-Val125
Accession	Q9GZQ8
Calculated Molecular Weight	14.8 kDa
Observed molecular weight	15 kDa
Tag	None

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 2mM DTT, pH 8.0 . Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Microtubule-associated proteins 1A/1B light chain 3B (MAP1LC3B) is a member of the highly conserved ATG8 protein family. ATG8 proteins are present in all known eukaryotic organisms. MAP1LC3B is one of the four genes in the MAP1LC3 subfamily (others include MAP1LC3A, MAP1LC3C, and MAP1LC3B2). It is most abundantly expressed in heart, brain, skeletal muscle and testis. LMAP1LC3B is a subunit of neuronal microtubule and functions in formation of

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human MAP1LC3B Protein

Catalog Number:PKSH033588



autophagosomal vacuoles (autophagosomes). It associated MAP1A and MAP1B proteins, which are involved in microtubule assembly and important for neurogenesis. MAP1LC3B also plays a role in autophagy, a process that involves the bulk degradation of cytoplasmic component.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com