

Recombinant Human JAM-A/F11R Protein (His Tag)

Catalog Number:PKSH033607



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Junctional Adhesion Molecule A;JAM-A;Junctional Adhesion Molecule 1;JAM-1;Platelet F11 Receptor;Platelet Adhesion Molecule 1;PAM-1;CD321;F11R;JAM1;JCAM;JAMA;JCAM;KAT
Species	Human
Expression Host	HEK293 Cells
Sequence	Ser28-Val238
Accession	Q9Y624
Calculated Molecular Weight	23.9 kDa
Observed molecular weight	28-31 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, 100mM Glycine, pH 7.5. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer informati
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Junctional Adhesion Molecule A (JAM-A) is a single-pass type I membrane protein that belongs to the immunoglobulin superfamily. JAM-A contains 2 Ig-like V-type (immunoglobulin-like) domains and Interacts with the ninth PDZ domain. JAM-A is localized to the tight junctions of both epithelial and endothelial cells. JAM-A seems to be involved in epithelial tight junction formation. JAM-A appears early in primordial forms of cell junctions and recruits PARD3. The association of the PARD6-PARD3 complex may prevent the interaction of PARD3 with JAM-A, thereby preventing tight junction

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human JAM-A/F11R Protein (His Tag)

Catalog Number:PKSH033607



assembly. JAM-A plays a role in regulating monocyte transmigration involved in regulating integrity of the epithelial barrier. In the case of orthoreovirus infection, JAM-A serves as receptor for the virus.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com