

Recombinant Human HMGB1 Protein

Catalog Number:PKSH033654



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	High Mobility Group Protein B1;High Mobility Group Protein 1;HMG-1;HMGB1;HMG1
Species	Human
Expression Host	E.coli
Sequence	Met1-Phe89
Accession	P09429
Calculated Molecular Weight	10.4 kDa
Observed molecular weight	14 kDa
Tag	None

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 50mM HEPES, 500mM NaCl, 0.5mMDTT, pH 7.9. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual for detailed information.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

High mobility group protein B1 is a member of the HMGB family consisting of three members; HMGB1; HMGB2 and HMGB3. It Contains 2 HMG box DNA-binding domains entitled box A and box B and It is a highly negative-charged C terminus. As a nuclear protein; HMGB1 stabilizes nucleosomes and allows bending of DNA that facilitates gene transcription which is essential for individual survival. Meanwhile; it is revealed that HMGB1 can also act as a cytokine extracellularly and regulates monocyte; T cell; dendritic cell activities in inflammatory responses.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017