

Recombinant Human SerpinE2/SERPINE2 Protein (His Tag)



Catalog Number:PKSH033671

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Glia-derived Nexin;GDN;Peptidase inhibitor 7;PI-7;Protease nexin 1;PN-1;Protease nexin I;Serpin E2;SERPINE2;PI7;PN1
Species	Human
Expression Host	HEK293 Cells
Sequence	Ser20-Pro397
Accession	P07093-2
Calculated Molecular Weight	42.7 kDa
Observed molecular weight	50 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < - 20°C.
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0.
Reconstitution	Not Applicable

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Serpin E2 is a secreted protein that belongs to the serpin family. Serpin E2 is a serine protease inhibitor with activity toward thrombin, trypsin, and urokinase. Serpin E2 expression is weak or absent in all normal pancreas and chronic pancreatitis tissue. In contrast, it was strongly over-expressed in the majority of pancreatic carcinoma as well as gastric and colorectal cancer samples. Serpin E2 promotes neurite extension by inhibiting thrombin. It also can bind heparin. It has been shown that Serpin E2 is a novel target of ERK signaling involved in human colorectal tumorigenesis. It plays an important role in controlling male fertility because its knockout male mice show a marked impairment in fertility from the onset of sexual maturity and its abnormal expression is found in the semen of men with seminal dysfunction.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017