

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	follistatin isoform FST317;Follistatin;FS;FSActivin-binding protein;FST
Species	Human
Expression Host	HEK293 Cells
Sequence	Gly30-Asn317
Accession	P19883
Calculated Molecular Weight	32.4 kDa
Observed molecular weight	33-42 kDa
Tag	C-His

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

Background

Follistatin 288 is a secreted glycoprotein that was first identified as a follicle-stimulating hormone inhibiting substance in ovarian follicular fluid. Human follistatin 288 cDNA encodes a 317 amino acid (aa) protein with a 29aa signal sequence, and a 288 aa mature region. Follistatin shows the highest affinity for activins due to its extended configuration. Genetic deletion of follistatin in mice, or expression of only the Follistatin form, is perinatally lethal due to defects of lung, skin and musculoskeletal system. Follistatins also regulate hematopoietic stem cell adhesion to fibronectin via FS2.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017