

Recombinant Human EphB1/EPHT2 Protein (aa 564-984, His Tag)



Catalog Number:PKSH033688

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

| | |
|------------------------------------|--|
| Synonyms | Ephrin Type-B Receptor 1;ELK;EPH Tyrosine Kinase 2;EPH-Like Kinase 6;EK6;hEK6;Neuronally-Expressed EPH;Related Tyrosine Kinase;NET;Tyrosine-Protein Kinase Receptor EPH-2;EPHB1;ELK;EPHT2;HEK6 |
| Species | Human |
| Expression Host | HEK293 Cells |
| Sequence | Ser564-Ala984 |
| Accession | AAI11745.1 |
| Calculated Molecular Weight | 48.8 kDa |
| Observed molecular weight | 35-50 kDa |
| Tag | C-His |
| Bioactivity | Immobilized Human EphB1-His at 10µg/ml(100 µl/well) can bind Mouse EFN2-Fc(Cat: PKSM041012). The ED50 of Human EphB1-His is 53. 1 ug/ml. |

Properties

| | |
|-----------------------|--|
| Purity | > 95 % as determined by reducing SDS-PAGE. |
| Endotoxin | < 1.0 EU per µg of the protein as determined by the LAL method. |
| Storage | Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months. |
| Shipping | This product is provided as lyophilized powder which is shipped with ice packs. |
| Formulation | Lyophilized from a 0.2 µm filtered solution of 20mM Tris-HCl, 150mM NaCl, pH 8.0. Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the print |
| Reconstitution | Please refer to the printed manual for detailed information. |

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Ephrin Type-B Receptor 1 (EPHB1) is a single-pass type I membrane protein that belongs to the Ephrin-B family of receptor tyrosine kinases that is involved in embryonic nervous and vascular system development. EPHB1/EPHT2 contains two fibronectin type-III domains; one protein kinase domain and one SAM (sterile α motif) domain. EPHB1

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human EphB1/EPHT2 Protein (aa 564-984, His Tag)



Catalog Number:PKSH033688

could stimulate fibroblast motility on extracellular matrix in a kinase-dependent manner; which also correlated with its association with Grb7; an adaptor molecule implicated in the regulation of cell migration. It binds to ephrin-B1; ephrin-B2 and ephrin-B3. EPHB1 plays an important roles in diverse biological processes including nervous system development; angiogenesis; and neural synapsis formation and maturation and may be involved in cell-cell interactions in the nervous system.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017