

Recombinant Human CFHR4 Protein (His Tag)

Catalog Number:PKSH033710



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Complement factor H-related protein 4;FHR-4;CFHR4;CFHL4
Species	Human
Expression Host	HEK293 Cells
Sequence	Glu20-Glu578
Accession	Q92496
Calculated Molecular Weight	64.1 kDa
Observed molecular weight	90-110 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Complement factor H-related protein 4 is a protein that in humans is encoded by the CFHR4 gene. It is a secreted protein that belongs to the complement factor H protein family. Members of the H-related protein family are exclusively composed of individually folded protein domains. CFHR4 is synthesized as a 578 amino acid precursor that contains an 19 amino acid signal peptide and a 569 amino acid mature chain. Human CFHR4 is involved in complement regulation. It can associate with lipoproteins and may play a role in lipid metabolism.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017