

Recombinant Mouse IL-21R Protein (His Tag)

Catalog Number:PKSM041316



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Interleukin-21 receptor;IL-21 receptor;IL-21R;Il21r;Nilr
Species	Mouse
Expression Host	HEK293 Cells
Sequence	Cys20-Pro236
Accession	Q9JHX3
Calculated Molecular Weight	25.7 kDa
Observed molecular weight	34-38 kDa
Tag	C-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

Interleukin-21 receptor (IL-21R) is a type I transmembrane glycoprotein within the class I cytokine receptor family, type 4 subfamily. IL-21R is expressed mainly on B cells, NK cells, and activated T cells, but is also found on dendritic cells, alternatively activated macrophages, intestinal lamina propria fibroblasts and epithelial cells, and keratinocytes. Both IL-21 and IL-4 are necessary for efficient B cell IgG1 production and normal germinal center architecture. B cell IL-21 R engagement induces Blimp-1(which mediates plasma cell differentiation), and is important for memory responses. IL-21R engagement on mouse NK cells enhances their cytotoxic activity and IFN- γ production. IL-21R engagement on CD8+ T cells aids control of viral infection and tumor growth? IL-21R is also necessary for sufficient numbers of regulatory T cells to combat chronic inflammation. IL-21R expression is often up-regulated in allergic skin inflammation, systemic

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Mouse IL-21R Protein (His Tag)

Catalog Number:PKSM041316



lupus erythematosus and diffuse large B cell lymphoma (DLBCL).

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com