

Recombinant Human ABHD4 Protein (His Tag)

Catalog Number:PKSH031257



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	ABH4
Species	Human
Expression Host	Baculovirus-Insect Cells
Sequence	Met 1-Asp 342
Accession	NP_071343.2
Calculated Molecular Weight	41.0 kDa
Observed molecular weight	40 kDa
Tag	N-His

Properties

Purity	> 82 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile 50mM Tris, 100mM NaCl, pH 8.0 Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 82 % as determined by reducing SDS-PAGE.

Background

Abhydrolase domain containing 4 (ABHD4), also known as alpha/beta-hydrolase 4 (ABH4), or lyso-N-acylphosphatidylethanolamine lipase, which belongs to the ABHD4/ABHD5 subfamily of peptidase S33 family. Abhydrolase domain containing (ABHD) gene was a small group belongs to alpha/beta hydrolase superfamily. Known members of this group are all found to be involved in important biochemical processes and related to various diseases. The alpha/beta-hydrolase 4 (ABH4) is a lysophospholipase/phospholipase B that selectively hydrolyzes N-acyl phosphatidylethanolamines (NAPEs) and lysoNAPEs. ABH4 accepts lysoNAPEs bearing both saturated and polyunsaturated N-acyl chains as substrates and displays a distribution that closely mirrors lysoNAPE-lipase activity in mouse tissues. The existence of an NAPE-PLD-independent route for NAE biosynthesis and suggest that ABH4 plays a

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Human ABHD4 Protein (His Tag)

Catalog Number:PKSH031257



role in this metabolic pathway by acting as a (lyso)NAPE-selective lipase.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com