

# Recombinant Human cIAP1/HiAP2 Protein (AVI Tag)

Catalog Number:PKSH031263



**Note:** Centrifuge before opening to ensure complete recovery of vial contents.


## Description

<b>Synonyms</b>	API1;c-IAP1;cIAP1;Hiap-2;HIAP2;MIHB;RNF48
<b>Species</b>	Human
<b>Expression Host</b>	E.coli
<b>Sequence</b>	Glu 144-Leu 356
<b>Accession</b>	NP_001157.1
<b>Calculated Molecular Weight</b>	26.5 kDa
<b>Observed molecular weight</b>	26.5 kDa
<b>Tag</b>	C-Avi
<b>Bioactivity</b>	Measured by its ability to inhibit DEVD-AFC cleavage activity in cell extracts activated by addition of cytochrome c and dATP. The IC50 for this effect is typically 25-750 nM.

## Properties

<b>Purity</b>	> 92 % as determined by reducing SDS-PAGE.
<b>Endotoxin</b>	Please contact us for more information.
<b>Storage</b>	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
<b>Shipping</b>	This product is provided as lyophilized powder which is shipped with ice packs.
<b>Formulation</b>	Lyophilized from sterile 10mM Tris, 5% glycerol, 0.5mM EDTA, 5mM DTT, pH 7.5 Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed ma
<b>Reconstitution</b>	Please refer to the printed manual for detailed information.

## Data



> 92 % as determined by reducing SDS-PAGE.

## Background

The cellular inhibitor of apoptosis protein-1 (cIAP1) is a member of the Inhibitor of Apoptosis family proteins (IAP) whose members are characterized by a novel domain of about 70 amino acids termed baculoviral IAP repeats (BIRs). The BIR domains of cIAP1 and cIAP2 bind to caspases, the key effector proteases of apoptosis. The IAP protein family which can enhance cell survival are crucial regulators of programmed cell death. Both cIAP1 and cIAP2 are the E3

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: [www.elabscience.com](http://www.elabscience.com)

Tel: 1-832-243-6086

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)

Fax: 1-832-243-6017

# Recombinant Human cIAP1/HiAP2 Protein (AVI Tag)

Catalog Number:PKSH031263



ubiquitin protein isopeptide ligases for Smac, taking part in promoting cancer survival through functioning as E3 ubiquitin ligases. Removal of cIAP1 by genetic deletion may result in NF- $\kappa$ B signaling activation that induces TNF $\alpha$  production and in killing sensitive tumor cells through enhanced TNF-R1 death-receptor signaling and caspase 8 activation. The substrate-dependent E3 activity of cIAPs is mediated by their RING domains and is dependent on the specific interactions between cIAPs and Smac. cIAP1 and cIAP2 are also reported to be regulators of NF- $\kappa$ B activation upon TNF treatment.

---

## For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: [www.elabscience.com](http://www.elabscience.com)

Email: [techsupport@elabscience.com](mailto:techsupport@elabscience.com)