

Recombinant Human ANXA5/Annexin V/Annexin A5 Protein

Catalog No. PKSH031601

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	ANXA5;Annexin V;Unconjugated Annexin V;Annexin-5;Anchorin CII;Calphobindin I;Endonexin II;Lipocortin V;Placental anticoagulant protein 4;Placental anticoagulant protein I;Thromboplastin inhibitor;Vascular anticoagulant-alpha
Species	Human
Expression Host	E.coli
Sequence	Met 1-Asp 320
Accession	P08758
Calculated Molecular Weight	36.8 kDa
Observed molecular weight	37 kDa
Tag	C-His
Bioactivity	Not validated for activity

Properties

Purity	> 98 % as determined by reducing SDS-PAGE.
Endotoxin	< 0.1 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 98 % as determined by reducing SDS-PAGE.

For Research Use Only

Background

Annexin V predicts a molecular mass of 36.8 kDa. is commonly used to detect apoptotic cells by its ability to bind to phosphatidylserine, a marker of apoptosis when it is on the outer leaflet of the plasma membrane. Annexin V forms a shield around negatively charged phospholipid molecules. The formation of an annexin V shield blocks the entry of phospholipids into coagulation reactions.