

Recombinant Mouse MCPT1 Protein (His Tag)

Catalog Number:PKSM040796



Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	AV080368;Mcp-1
Species	Mouse
Expression Host	HEK293 Cells
Sequence	Met 1-Lys 246
Accession	NP_032596.1
Calculated Molecular Weight	26.8 kDa
Observed molecular weight	32-34 kDa
Tag	C-His

Properties

Purity	> 97 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4 Normally 5 % - 8 % trehalose, mannitol and 0.01 % Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 97 % as determined by reducing SDS-PAGE.

Background

Mast Cell Protease 1 (MMCP-1), also known as MCP-1, MCPT-1 and β -chymase, is a member of the Chymase family of chymotrypsin-like serine proteases. MCPT-1 is a 26 kDa β -chymase that is a component of mast cell granules. It is a 226 amino acid (aa) protein that has a conserved pattern of six cysteines and one potential glycosylation site. The granule-derived mouse mast cell proteases-1 and -2 (mMCP-1 and -2) colocalize in similar quantities in mucosal mast cells but micrograms of mMCP-1 compared with nanograms of mMCP-2 are detected in peripheral blood during intestinal nematode infection. mMCP-1 isolated from serum is complexed with serpins and both the accumulation and the longevity of mMCP-1 in blood is due to complex formation, protecting it from a pathway that rapidly clears mMCP-2, which is unable to form complexes with serpins. Expression of mMCP-1 is largely restricted to intraepithelial MMC and is thought

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017

Recombinant Mouse MCPT1 Protein (His Tag)

Catalog Number:PKSM040796



to play a role in the regulation of epithelial permeability. Its activation is completed by the removal of a two residue N-terminal propeptide by a dipeptidyl peptidase (Cathepsin C). MCPT-1 is upregulated in the intestine in response to nematode infection, or in systemic mucosa in response to anaphylaxis. Like human α -chymase, MCPT-1 is capable of the conversion of angiotensin I to angiotensin II, which plays a key role in the regulation of arterial pressure. The intestinal inflammation associated with gastrointestinal helminths is partly mediated by mMCP-1.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Tel: 1-832-243-6086

Fax: 1-832-243-6017

Web: www.elabscience.com

Email: techsupport@elabscience.com