

Recombinant DENV (type 4, strain Philippines/H241/1956) E / Envelope Protein (ECD, His Tag)

Catalog No. PKSV030130

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	E Protein;DENV;Dengue virus;type 4;strain Philippines/H241/1956
Species	DENV
Expression Host	Baculovirus-Insect Cells
Sequence	Met280-Gly674
Accession	AAX48017.1
Calculated Molecular Weight	44.8 kDa
Tag	C-His

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile 20 mM Tris, 500 mM NaCl, 10% glycerol, pH 8.0. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

For Research Use Only