

Recombinant DENV-2 (Strain New Guinea C) Capsid protein / DENV-C Protein (His Tag)

Catalog No. PKSV030132

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Capsid protein;DENV-C Protein;Strain New Guinea C;Dengue virus
Species	DENV
Expression Host	E.coli
Sequence	Met1-Arg100
Accession	AHF52833.1
Calculated Molecular Weight	14 kDa
Tag	N-His

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	Please contact us for more information.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

Background

death-domain associated protein

For Research Use Only