

Recombinant EBOV (subtype Bundibugyo, strain Uganda 2007) VP40 / Matrix protein VP40 Protein (His Tag)

Catalog No. PKSV030139

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	EBOV VP40;MatrixVP40
Species	EBOV
Expression Host	E.coli
Sequence	Ala31-Lys326
Accession	YP_003815434.1
Calculated Molecular Weight	34.4 kDa.
Tag	N-His

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	Please contact us for more information.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile 50 mM Tris, 50 mM NaCl, pH 8.0. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

For Research Use Only