

Recombinant Hantaan virus HTNV (strain 84FLi) Nucleocapsid / NP Protein (His Tag)

Catalog No. PKSV030173

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	NP Protein;HTNV
Species	HTNV
Expression Host	Baculovirus-Insect Cells
Sequence	Met1-Leu429
Accession	AAK66740.1
Calculated Molecular Weight	49.7 kDa.
Tag	C-His

Properties

Purity	> 80 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Store at < -20°C, stable for 6 months. Please minimize freeze-thaw cycles.
Shipping	This product is provided as liquid. It is shipped at frozen temperature with blue ice/gel packs. Upon receipt, store it immediately at < -20°C.
Formulation	Supplied as sterile 20 mM Tris, 500 mM NaCl, 20 % glycerol, pH 8.0.
Reconstitution	Not Applicable

Data



> 80 % as determined by reducing SDS-PAGE.