

Recombinant ZIKV (strain Zika SPH2015) NS1 protein (C-His Tag)



Catalog Number:PKSV030273

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	Zika virus;NS1 protein;strain Zika SPH2015
Species	ZIKV
Expression Host	Baculovirus-Insect Cells
Sequence	Val796-Leu1157
Accession	ALU33341.1
Calculated Molecular Weight	42.6 kDa
Tag	C-His

Properties

Purity	> 90 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile 20 mM Tris, 500 mM NaCl, pH 8.0, 10% glycerol. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 90 % as determined by reducing SDS-PAGE.

Background

Zika virus NS1 antigen is one of seven non-structural proteins. NS1 is involved in RNA replication. The possible effects of NS1 on hosts include: localization to host cell surface and secreted extracellularly, modulates signalling of innate immune system, has possible damages to platelets and endothelial cells through anti-NS1 antibodies.

For Research Use Only

A Reliable Research Partner in Life Science and Medicine

Toll-free: 1-888-852-8623

Web: www.elabscience.com

Tel: 1-832-243-6086

Email: techsupport@elabscience.com

Fax: 1-832-243-6017