

Recombinant Human G-CSFR/CD114/CSF3R (C-Fc)

Catalog No. PKSH033845

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Synonyms	CD114 antigen;CD114;colony stimulating factor 3 receptor (granulocyte);CSF3R;Csfgr;G-CSF R;G-CSF receptor;GCSFR;G-CSFR;GCSFRG-CSF-R
Species	Human
Expression Host	HEK293 Cells
Sequence	Glu25-His627
Accession	Q99062
Calculated Molecular Weight	93.7 kDa
Observed molecular weight	110-130 kDa
Tag	C-Fc

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Background

Granulocyte Colony Stimulating Factor Receptor (G-CSFR), also known as CD114, the protein encoded by this gene is the receptor for colony stimulating factor 3, a cytokine that controls the production, differentiation, and function of granulocytes. The encoded protein, which is a member of the family of cytokine receptors, may also function in some cell surface adhesion or recognition processes. Mutations in the G-CSF receptor leading to carboxy-terminal truncation transduce hyperproliferative growth responses, and are implicated in the pathological progression of severe congenital neutropenia (SCN) to acute myelogenous leukemia (AML). Additionally, autocrine/paracrine stimulation of G-CSFR may be important in the biology of solid tumors, including metastasis.

For Research Use Only