

Recombinant Human HLA-A*0201 GP100 complex Protein (C-Avi-10His) Biotinylated

Catalog No. PKSH034058

Note: Centrifuge before opening to ensure complete recovery of vial contents.


Description

Synonyms	B2M;HLA-A
Species	Human
Expression Host	HEK293 Cells
Sequence	Ile21-Met119&Gly25-Ile308
Accession	P61769&P01892
Calculated Molecular Weight	50.9 kDa
Observed molecular weight	55-60 kDa
Tag	C-Avi-His

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	< 1.0 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 500mM NaCl, pH 7.4. Normally 5% - 8% trehalose, mannitol and 0.01% Tween 80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Data



> 95 % as determined by reducing SDS-PAGE.

Background

GP100 Plays a central role in the biogenesis of melanosomes. Involved in the maturation of melanosomes from stage I to II. The transition from stage I melanosomes to stage II melanosomes involves an elongation of the vesicle, and the

For Research Use Only

appearance within of distinct fibrillar structures.