

## VERO C1008 [Vero E6] Cells Complete Medium

Cat. No: CM-0491  
Size: 125mL×4

### General Information

Product Form	Liquid
Concentration	Ready-to-use
Component	MEM( <a href="#">PM150467</a> )+Nutrients+1% P/S( <a href="#">PB180120</a> )
Bacterial detection	Negative
Fungal detection	Negative
Mycoplasma detection	Negative
Endotoxin content	≤3
Storage	2-8°C ,Shading Light
Shipping	Ice bag
Expiration date	3 months

### Background

VERO C1008 [Vero E6] Cell Complete Medium is a ready-to-use medium that can be directly used for the culture of VERO C1008 [Vero E6]cell

### Notes

1. This product is for research use only.;
2. This product is sterilized by 0.1μm filtration.;
3. It is necessary to pay attention to the aseptic operation and avoid the contamination during the culture.;
4. It is not suitable for long time storage at room temperature.;
5. This product is a ready-to-use medium. If there is no special need, don't add serum, penicillin and streptomycin. It can be used directly.