

Recombinant Human G-CSF protein(His Tag)

Catalog No. PKSH034164

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Synonyms	CSF-3;MGI-1G;GM-CSF beta;pluripoietin
Species	Human
Expression Host	E.coli
Sequence	Thr 31-Pro 204
Accession	NP_757373
Calculated Molecular Weight	19.5 kDa
Observed molecular weight	21 kDa
Tag	N-His

Properties

Purity	> 98 % as determined by reducing SDS-PAGE.
Endotoxin	< 0.1 EU per µg of the protein as determined by the LAL method.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	Please refer to the printed manual for detailed information.

Background

Granulocyte/macrophage colony-stimulating factors are cytokines that act in hematopoiesis by controlling the production, differentiation, and function of 2 related white cell populations of the blood, the granulocytes and the monocytes-macrophages. This CSF induces granulocytes.

For Research Use Only