

Recombinant Human Biglycan/BGN protein (His tag)

Catalog No. PDMH100070

Note: Centrifuge before opening to ensure complete recovery of vial contents.

Description

Synonyms	Biglycan;Bone/Cartilage Proteoglycan I;PG-S1;BGN;SLRR1A
Species	Human
Expression Host	HEK293 Cells
Sequence	Glu20-Lys368
Accession	P21810
Calculated Molecular Weight	40.4 kDa
Observed molecular weight	42 kDa
Tag	C-His
Bioactivity	Not validated for activity

Properties

Purity	> 95 % as determined by reducing SDS-PAGE.
Endotoxin	Please contact us for more information.
Storage	Generally, lyophilized proteins are stable for up to 12 months when stored at -20 to -80°C. Reconstituted protein solution can be stored at 4-8°C for 2-7 days. Aliquots of reconstituted samples are stable at < -20°C for 3 months.
Shipping	This product is provided as lyophilized powder which is shipped with ice packs.
Formulation	Lyophilized from sterile PBS, pH 7.4. Normally 5 % - 8 % trehalose, mannitol and 0.01% Tween80 are added as protectants before lyophilization. Please refer to the specific buffer information in the printed manual.
Reconstitution	It is recommended that sterile water be added to the vial to prepare a stock solution of 0.5 mg/mL. Concentration is measured by UV-Vis

Background

Biglycan is a 200-350 kD proteoglycan consisting of a 45 kD core protein and two chondroitin/dermatan sulfate glycosaminoglycan chains. Biglycan binds to TGF- β . It also binds to collagen type I in low ionic strength (less than 3 mM phosphate) buffer. At higher ionic strengths, Biglycan does not bind to collagen type I. It enhances the inhibition effect of TGF- β on osteoclast proliferation at a concentration of 4-20 mg/ml. It also prevents the attachment of CHO cells to fibronectin, with a 50% inhibition at 17-21 mg/ml.

For Research Use Only